

Examination Couch €100

VERSATILE • COMFORTABLE • CONVENIENT



CONVENIENCE for the care-provider
COMFORT for the Patient

Examination Couch



Electrically Operated
Trendelenburg
5" (12.7cm) Santoprene Total Lock Castors
(2 Brake, 2 Swivel)
Mattress density 50mm

Technical Specification

Minimum Height 500 mm
Maximum Height 1110mm
Total length 1180mm
Total width 570mm
Trendlenberg upto 45° (headside only)



The EC100 is available in a choice of mattress and frame colours to match your labour room colour scheme.
If your practice has any special requirement, call the R&D department at **PHOENIX**. They will be happy to provide you a solution



CDSO Registration No. : 000387 / 2020

Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Tamil Nadu – 600132, India

Toll Free India : 1800 425 3503, Phone: +91 44 2478 1705 / 2227,

✉ India : sales@pmsind.com, International: exports@pmsind.com

🌐 www.phoenixmedicalsystems.com