

Electronic Baby Weighing Scale BWS101



SPECIFICATIONS

Accuracy.....5 g

Resolution.....1 g

Weight

Maximum weight.....10 kg

Minimum weight10 g

Display

Display type.....LCD display

Temperature

Operating Temperature.....-10°C to 50°C

Power requirements

Input.....230 V AC, 50 Hz

Battery backup.....6 V, 1.2 Ah

Battery working time.....1 hour (when fully charged)

Battery charging time.....2 hours (approximate charged)

Dimensions of base

Length.....300 mm (12 inches)

Width265 mm (10 inches)

Height.....85 mm (3.5 inches)

Pan

Length.....510 mm (20 inches)

Width300 mm (12 inches)

Height.....85 mm (3.5 inches)

Material.....Transparent acrylic

Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Tamil Nadu – 600132, India

Toll Free India : 1800 425 3503, Phone: +91 44 2478 1705 / 2227,

✉ India : sales@pmsind.com, International: exports@pmsind.com

 www.phoenixmedicalsyste.ms.com