

Face Shield Specifications

FACE SHIELDS TO PROTECT HEALTH WORKERS FROM CORONAVIRUS

Salient Features

- Clear & Transparent Front Shield
- Quick & Easy To Assemble
- Environmentally Friendly Sponge
- Can Be Disinfected & Reused
- Adjustable Strap For Convenience
- Affordable Price – INR 50.00 / Shield



Technical Specifications

PARAMETER	DESCRIPTION
Shield Material	PET-G
Foam Material	Polyurethane
Elastic Strap	Polyester Fibre
Disinfection	With Isopropanol (70%)
Length	350 mm x 15 mm
Shield Size	360 mm x 180 mm
Temperature	< 55 degrees C
Transparency	100% Clear