

Specifications For Infra Red Thermometer



Forehead Non-Contact IR Thermometer

Model Name	AshaFit™
Model No	AF-XX20
Response Time	2 seconds
Measuring Distance	3-5cm
Body Temperature Measurement Range	35°C - 42°C / 95°F-107°F
Power Source	2 AAA Batteries
Accuracy	±0.2°C / ±0.5°F
Resolution ratio	±0.1 °C / °F
Dimensions	150mm x 100mm x 43mm
Display	OLED Display
°C / °F Selectable	Yes
Low Battery Indication	Yes
Alarm Sound	Yes
Operating Temperature	16°C-40°C
Humidity	<85% RH
Warranty	1 Year