



PERSONAL PROTECTIVE EQUIPMENT KIT

QUALITY STANDARDS

Coverall : SITRA, ASTM F 1670, ISO 16603

3-Ply Mask : SITRA, ASTM F 1670

Hand Gloves: CE Certified

Goggles: Comply to EN166:2001 standard



SITRA
Service for Excellence

The Personal Protective Equipment kit is designed to protect you from harmful viruses, bacteria and other infectious agents or fluids. The PPE kit contains all essentials to avoid any bodily exposure to such infections.