

ASHAPuls™



₹2500/-

Portable Pulse Oximeter

Model Name	AshaPuls™
Model No	AO-XX20
Measurement Technology	Reflective
SpO2 measuring range	80-99%
SpO2 accuracy	80-99%, ±3%
Pulse rate measuring range	40-200 bpm
Pulse rate accuracy	a. 40-100bpm: ±5bpm b. 100-200bpm: ±5 %
Power Source	2 AAA Batteries
Resolution ratio	SpO2: 1 %, Pulse : 1 bpm
Low Battery Indication	Yes
Operating Temperature	16°C-40°C
Humidity	<85% RH
Warranty	1 Year

इयावइथ्या लडोन डरे कवडल वडुड ASHA PULS™ के इवाथे

कंरुंनूडुडुनं रिसरुसु लेडुडुसु डु.लु.

GreenOcean Research Labs Pvt Ltd

#76, 7th Street, Poompuhar Nagar, Kolathur, Chennai - 600 099; Ph: +91 44 2650 0137

Factory: MU-E, 'Avyakthananda', AMTZ campus, Pragati Maidan, Visakhapatnam- 5300031