

ASHAPuls™



₹3500/-

Wireless Pulse Oximeter

Model Name	AshaPuls™
Model No	AO-XX20W
Measurement Technology	Reflective
SpO2 measuring range	80-99%
SpO2 accuracy	80-99%, ±2%
Pulse rate measuring range	40-200 bpm
Pulse rate accuracy	a. 40-100bpm: ±3bpm b. 100-200bpm: ± 3 %
Wireless communication	Bluetooth Low Energy
Device Compatibility	iOS, Android
Power Source	Rechargeable Battery, 3.7V Li-Po, 250mAh
Resolution ratio	SpO2: 1 %, Pulse : 1 bpm
Low Battery Indication	Yes
Operating Temperature	16°C-40°C
Humidity	<85% RH
Warranty	1 Year

इयावइथिया लडोन डेर कवडाल ववडुड ASHAPuls™ के इवाथे

कंरुंनुवडुडन रिसरुसु लेडुडनु डी.ली.

GreenOcean Research Labs Pvt Ltd

#76, 7th Street, Poompuhar Nagar, Kolathur, Chennai - 600 099; Ph: +91 44 2650 0137

Factory: MU-E, 'Avyakthananda", AMTZ campus, Pragati Maidan, Visakhapatnam- 5300031